



CL070 rabbit pAb

Cat No.:ES19439

For research use only

Overview

Product Name	CL070 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CL070 AA range: 232-282
Specificity	This antibody detects endogenous levels of human CL070
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Uncharacterized protein C12orf70
Gene Name	C12orf70
Cellular localization	Membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	341346
Human Swiss-Prot Number	A6NFE2
Alternative Names	
Background	

