

# NCKPL rabbit pAb

Cat No.:ES19441

For research use only

## Overview

Product Name	NCKPL rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NCKPL AA range: 582-632
Specificity	This antibody detects endogenous levels of human NCKPL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Nck-associated protein 1-like (Hematopoietic protein 1) (Membrane-associated protein HEM-1)
Gene Name	NCKAP1L HEM1
Cellular localization	Cell membrane ; Single-pass membrane protein ; Cytoplasmic side . Cytoplasm . Localizes to the leading edge of polarized neutrophils. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3071
Human Swiss-Prot Number	P55160
Alternative Names	
Background	

