



UBR4 rabbit pAb

Cat No.:ES19443

For research use only

Overview

Product Name	UBR4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human UBR4 AA range: 809-859
Specificity	This antibody detects endogenous levels of human UBR4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	E3 ubiquitin-protein ligase UBR4 (EC 6.3.2.-) (600 kDa retinoblastoma protein-associated factor) (N-recognin-4) (Retinoblastoma-associated factor of 600 kDa) (RBAF600) (p600) (Zinc finger UBR1-type pr
Gene Name	UBR4 KIAA0462 KIAA1307 RBAF600 ZUBR1
Cellular localization	Membrane ; Multi-pass membrane protein . Cytoplasm. Cytoplasm, cytoskeleton. Nucleus. Concentrates at the leading edge of membrane structures involved in actin motility.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23352
Human Swiss-Prot Number	Q5T4S7
Alternative Names	
Background	

