



CL18C rabbit pAb

Cat No.:ES19444

For research use only

Overview

Product Name	CL18C rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CL18C AA range: 272-322
Specificity	This antibody detects endogenous levels of human CL18C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	C-type lectin domain family 18 member C (Mannose receptor-like protein 3)
Gene Name	CLEC18C MRLP3
Cellular localization	Secreted . Endoplasmic reticulum . Golgi apparatus . Endosome .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	283971
Human Swiss-Prot Number	Q8NCF0
Alternative Names	
Background	

