



COFA1 rabbit pAb

Cat No.:ES19446

For research use only

Overview

Product Name	COFA1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human COFA1 AA range: 1231-1280
Specificity	This antibody detects endogenous levels of human COFA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Collagen alpha-1(XV) chain [Cleaved into: Endostatin (Endostatin-XV) (Related to endostatin) (Restin)]
Gene Name	COL15A1
Cellular localization	Secreted, extracellular space, extracellular matrix .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	1306
Human Swiss-Prot Number	P39059
Alternative Names	
Background	

