



CPEB4 rabbit pAb

Cat No.:ES19450

For research use only

Overview

Product Name	CPEB4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CPEB4 AA range: 667-717
Specificity	This antibody detects endogenous levels of human CPEB4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Cytoplasmic polyadenylation element-binding protein 4 (CPE-BP4) (CPE-binding protein 4) (hCPEB-4)
Gene Name	CPEB4 KIAA1673
Cellular localization	Cytoplasm . Cell projection, dendrite . Cell projection, dendritic spine . Cell junction, synapse, postsynaptic density . Cell projection, axon . Cell projection, growth cone . Endoplasmic reticulum . Cytoplasm, perinuclear region .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	80315
Human Swiss-Prot Number	Q17RY0
Alternative Names	
Background	

