



K1C27 rabbit pAb

Cat No.:ES19453

For research use only

Overview

Product Name	K1C27 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human K1C27 AA range: 18-68
Specificity	This antibody detects endogenous levels of human K1C27
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Keratin, type I cytoskeletal 27 (Cytokeratin-27) (CK-27) (Keratin-25C) (K25C) (Keratin-27) (K27) (Type I inner root sheath-specific keratin-K25irs3)
Gene Name	KRT27 KRT25C
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	342574
Human Swiss-Prot Number	Q7Z3Y8
Alternative Names	
Background	

