

EST1A rabbit pAb

Cat No.:ES19460

For research use only

Overview

Product Name	EST1A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EST1A AA range: 688-738
Specificity	This antibody detects endogenous levels of human EST1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Telomerase-binding protein EST1A (EC 3.1.-.-) (EST1-like protein A) (Ever shorter telomeres 1A) (Smg-6 homolog) (Telomerase subunit EST1A) (hSmg5/7a)
Gene Name	SMG6 C17orf31 EST1A KIAA0732
Cellular localization	Nucleus, nucleolus . Chromosome, telomere . Cytoplasm, cytosol . Particularly enriched in the nucleolus. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23293
Human Swiss-Prot Number	Q86US8
Alternative Names	
Background	

