

LUZP6 rabbit pAb

Cat No.:ES19464

For research use only

Overview

Product Name	LUZP6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LUZP6 AA range: 8-58
Specificity	This antibody detects endogenous levels of human LUZP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Leucine zipper protein 6 (Myeloproliferative disease-associated 6 kDa antigen)
Gene Name	LUZP6 MPD6
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	767558
Human Swiss-Prot Number	Q538Z0
Alternative Names	
Background	

