



PPDPF rabbit pAb

Cat No.:ES19470

For research use only

Overview

Product Name	PPDPF rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PPDPF AA range: 59-109
Specificity	This antibody detects endogenous levels of human PPDPF
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Pancreatic progenitor cell differentiation and proliferation factor (Exocrine differentiation and proliferation factor)
Gene Name	PPDPF C20orf149 EXPDF
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79144
Human Swiss-Prot Number	Q9H3Y8
Alternative Names	
Background	

