



UN13D rabbit pAb

Cat No.:ES19473

For research use only

Overview

| | |
|---------------------------------|--|
| Product Name | UN13D rabbit pAb |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human;Mouse;Rat |
| Recommended dilutions | WB 1:1000-2000 |
| Immunogen | Synthesized peptide derived from human UN13D AA range: 213-263 |
| Specificity | This antibody detects endogenous levels of human UN13D |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Protein Name | Protein unc-13 homolog D (Munc13-4) |
| Gene Name | UNC13D |
| Cellular localization | Cytoplasm. Membrane; Peripheral membrane protein. Late endosome. Recycling endosome. Lysosome. Colocalizes with cytotoxic granules at the plasma membrane. Localizes to endosomal exocytic vesicles. |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | 201294 |
| Human Swiss-Prot Number | Q70J99 |
| Alternative Names | |
| Background | |

