

# TM102 rabbit pAb

Cat No.:ES19474

For research use only

## Overview

Product Name	TM102 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TM102 AA range: 275-325
Specificity	This antibody detects endogenous levels of human TM102
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Transmembrane protein 102 (Common beta-chain associated protein) (CBAP)
Gene Name	TMEM102
Cellular localization	Cell membrane ; Single-pass membrane protein . Although it is found at the cell surface, the majority of the molecules seems to be located in intracellular compartments.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	284114
Human Swiss-Prot Number	Q8N9M5
Alternative Names	
Background	

