

# NAT6 rabbit pAb

Cat No.:ES19478

For research use only

## Overview

Product Name	NAT6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NAT6 AA range: 107-157
Specificity	This antibody detects endogenous levels of human NAT6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	N-acetyltransferase 6 (EC 2.3.1.-) (Protein fusion-2) (Protein fus-2)
Gene Name	NAT6 FUS2
Cellular localization	Cytoplasm, cytosol .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	24142
Human Swiss-Prot Number	Q93015
Alternative Names	
Background	

