

KHDR2 rabbit pAb

Cat No.:ES19481

For research use only

Overview

Product Name	KHDR2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KHDR2 AA range: 269-319
Specificity	This antibody detects endogenous levels of human KHDR2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	KH domain-containing, RNA-binding, signal transduction-associated protein 2 (Sam68-like mammalian protein 1) (SLM-1) (hSLM-1)
Gene Name	KHDRBS2 SLM1
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	202559
Human Swiss-Prot Number	Q5VWX1
Alternative Names	
Background	

