



HEMH rabbit pAb

Cat No.:ES19483

For research use only

Overview

Product Name	HEMH rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HEMH AA range: 169-219
Specificity	This antibody detects endogenous levels of human HEMH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Ferrochelatase, mitochondrial (EC 4.99.1.1) (Heme synthase) (Protoheme ferro-lyase)
Gene Name	FECH
Cellular localization	Mitochondrion inner membrane; Peripheral membrane protein; Matrix side.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	2235
Human Swiss-Prot Number	P22830
Alternative Names	
Background	

