



CECR1 rabbit pAb

Cat No.:ES19484

For research use only

Overview

Product Name	CECR1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CECR1 AA range: 243-293
Specificity	This antibody detects endogenous levels of human CECR1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Adenosine deaminase CECR1 (EC 3.5.4.4) (Cat eye syndrome critical region protein 1)
Gene Name	CECR1 ADA2 ADGF IDGFL
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	51816
Human Swiss-Prot Number	Q9NZK5
Alternative Names	
Background	

