



KAZD1 rabbit pAb

Cat No.:ES19486

For research use only

Overview

Product Name	KAZD1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KAZD1 AA range: 4-54
Specificity	This antibody detects endogenous levels of human KAZD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Kazal-type serine protease inhibitor domain-containing protein 1
Gene Name	KAZALD1 FKSG28 FKSG40 UNQ2945/PRO21184
Cellular localization	Secreted, extracellular space, extracellular matrix .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	81621
Human Swiss-Prot Number	Q96I82
Alternative Names	
Background	

