



WFKN1 rabbit pAb

Cat No.:ES19492

For research use only

Overview

Product Name	WFKN1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human WFKN1 AA range: 83-133
Specificity	This antibody detects endogenous levels of human WFKN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	WAP, kazal, immunoglobulin, kunitz and NTR domain-containing protein 1 (Growth and differentiation factor-associated serum protein 2) (GASP-2) (hGASP-2) (WAP, follistatin, immunoglobulin, kunitz and N
Gene Name	WFIKKN1 C16orf12 GASP2 WFIKKN PSEC0040
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	117166
Human Swiss-Prot Number	Q96NZ8
Alternative Names	
Background	

