

CL045 rabbit pAb

Cat No.:ES19493

For research use only

Overview

Product Name	CL045 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CL045 AA range: 45-95
Specificity	This antibody detects endogenous levels of human CL045
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Uncharacterized protein C12orf45
Gene Name	C12orf45
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	121053
Human Swiss-Prot Number	Q8N5I9
Alternative Names	
Background	

