

ARP3C rabbit pAb

Cat No.: ES19499

For research use only

Overview

Product Name ARP3C rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human ARP3C AA

range: 86-136

Specificity This antibody detects endogenous levels of human

ARP3C

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Actin-related protein 3C (Actin-related protein 11)

Gene Name ACTR3C ARP11

Cellular localization

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 653857 Human Swiss-Prot Number Q9C0K3

Alternative Names

+86-27-59760950

Background



23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.