

ME15P rabbit pAb

Cat No.:ES19500

For research use only

Overview

| | |
|--------------------------|--|
| Product Name | ME15P rabbit pAb |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human;Rat;Mouse; |
| Recommended dilutions | WB 1:1000-2000 |
| Immunogen | Synthesized peptide derived from human ME15P AA range: 66-116 |
| Specificity | This antibody detects endogenous levels of human ME15P |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20℃. Avoid repeated freeze-thaw cycles. |
| Protein Name | Putative methyltransferase-like protein 15P1 (EC 2.1.1.-) (Methyltransferase 5 domain-containing protein 2) (Methyltransferase-like protein 15 pseudogene 1) |
| Gene Name | METTL15P1 METT5D2 |
| Cellular localization | |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | |
| Human Swiss-Prot Number | POC7V9 |
| Alternative Names | |
| Background | |

