



TSPY4 rabbit pAb

Cat No.:ES19511

For research use only

Overview

Product Name	TSPY4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TSPY4 AA range: 90-140
Specificity	This antibody detects endogenous levels of human TSPY4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Testis-specific Y-encoded protein 4
Gene Name	TSPY4
Cellular localization	Cytoplasm . Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	728395
Human Swiss-Prot Number	POCV99
Alternative Names	
Background	

