



# PLPL3 rabbit pAb

Cat No.:ES19516

For research use only

## Overview

Product Name	PLPL3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PLPL3 AA range: 123-173
Specificity	This antibody detects endogenous levels of human PLPL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Patatin-like phospholipase domain-containing protein 3 (EC 3.1.1.3) (Acylglycerol O-acyltransferase) (EC 2.3.1.-) (Adiponutrin) (Calcium-independent phospholipase A2-epsilon) (iPLA2-epsilon)
Gene Name	PNPLA3 ADPN C22orf20
Cellular localization	Membrane ; Single-pass type II membrane protein . Lipid droplet .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	80339
Human Swiss-Prot Number	Q9NST1
Alternative Names	
Background	

