

SMRP1 rabbit pAb

Cat No.:ES19517

For research use only

Overview

Product Name	SMRP1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SMRP1 AA range: 130-180
Specificity	This antibody detects endogenous levels of human SMRP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Spermatid-specific manchette-related protein 1 (Ciliated bronchial epithelial protein 1) (Testis development protein NYD-SP22)
Gene Name	SMRP1 C9orf24 CBE1
Cellular localization	Cytoplasm, cytoskeleton . Nucleus . Cytoplasm . Localizes to the manchette in elongating spermatids (By similarity). Predominantly perinuclear in bronchial epithelial cells but also detected in the nucleus in some primary epithelial cells and in a number of cell lines. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84688
Human Swiss-Prot Number	Q8NCR6
Alternative Names	
Background	

