

# YAED1 rabbit pAb

Cat No.:ES19519

For research use only

## Overview

Product Name	YAED1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human YAED1 AA range: 18-68
Specificity	This antibody detects endogenous levels of human YAED1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Yae1 domain-containing protein 1
Gene Name	YAE1D1 C7orf36 GK003
Cellular localization	Cytoplasm . Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	57002
Human Swiss-Prot Number	Q9NRH1
Alternative Names	
Background	

