



MZT2A rabbit pAb

Cat No.:ES19523

For research use only

Overview

Product Name	MZT2A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MZT2A AA range: 72-122
Specificity	This antibody detects endogenous levels of human MZT2A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Mitotic-spindle organizing protein 2A (Mitotic-spindle organizing protein associated with a ring of gamma-tubulin 2A)
Gene Name	MZT2A FAM128A MOZART2A
Cellular localization	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Cytoplasm, cytoskeleton, spindle.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	653784
Human Swiss-Prot Number	Q6P582
Alternative Names	
Background	

