

# PDPK2 rabbit pAb

Cat No.:ES19524

For research use only

## Overview

Product Name	PDPK2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PDPK2 AA range: 254-304
Specificity	This antibody detects endogenous levels of human PDPK2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Putative 3-phosphoinositide-dependent protein kinase 2 (EC 2.7.11.1)
Gene Name	PDPK2
Cellular localization	Cytoplasm . Membrane ; Peripheral membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q6A1A2
Alternative Names	
Background	

