



C170L rabbit pAb

Cat No.:ES19526

For research use only

Overview

Product Name	C170L rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human C170L AA range: 13-63
Specificity	This antibody detects endogenous levels of human C170L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Cep170-like protein (CEP170 pseudogene 1)
Gene Name	CEP170P1 CEP170L KIAA0470L
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q96L14
Alternative Names	
Background	

