

NMNA2 rabbit pAb

Cat No.:ES19527

For research use only

Overview

Product Name	NMNA2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NMNA2 AA range: 218-268
Specificity	This antibody detects endogenous levels of human NMNA2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Nicotinamide mononucleotide adenylyltransferase 2 (NMN adenylyltransferase 2) (EC 2.7.7.1) (Nicotinate-nucleotide adenylyltransferase 2) (NaMN adenylyltransferase 2) (EC 2.7.7.18)
Gene Name	NMNAT2 C1orf15 KIAA0479
Cellular localization	Golgi apparatus membrane ; Lipid-anchor . Cytoplasmic vesicle membrane ; Lipid-anchor . Cytoplasm . Cell projection, axon . Delivered to axons with Golgi-derived cytoplasmic vesicles. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23057
Human Swiss-Prot Number	Q9BZQ4
Alternative Names	
Background	

