

# PHYIP rabbit pAb

Cat No.:ES19528

For research use only

## Overview

Product Name	PHYIP rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PHYIP AA range: 1-51
Specificity	This antibody detects endogenous levels of human PHYIP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Phytanoyl-CoA hydroxylase-interacting protein (Phytanoyl-CoA hydroxylase-associated protein 1) (PAHX-AP1) (PAHXAP1)
Gene Name	PHYHIP DYRK1AP3 KIAA0273
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9796
Human Swiss-Prot Number	Q92561
Alternative Names	
Background	

