



KIFA3 rabbit pAb

Cat No.:ES19529

For research use only

Overview

Product Name	KIFA3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KIFA3 AA range: 65-115
Specificity	This antibody detects endogenous levels of human KIFA3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Kinesin-associated protein 3 (KAP-3) (KAP3) (Smg GDS-associated protein)
Gene Name	KIFAP3 KIF3AP SMAP
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	22920
Human Swiss-Prot Number	Q92845
Alternative Names	
Background	

