



MBD5 rabbit pAb

Cat No.:ES19530

For research use only

Overview

Product Name	MBD5 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MBD5 AA range: 1057-1107
Specificity	This antibody detects endogenous levels of human MBD5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Methyl-CpG-binding domain protein 5 (Methyl-CpG-binding protein MBD5)
Gene Name	MBD5 KIAA1461
Cellular localization	[Isoform 1]: Nucleus. Chromosome. Associated with pericentric heterochromatin.; [Isoform 2]: Nucleus. Not associated with pericentric heterochromatin.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55777
Human Swiss-Prot Number	Q9P267
Alternative Names	
Background	

