

LIPJ rabbit pAb

Cat No.:ES19532

For research use only

Overview

| | |
|--------------------------|---|
| Product Name | LIPJ rabbit pAb |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human;Rat;Mouse; |
| Recommended dilutions | WB 1:1000-2000 |
| Immunogen | Synthesized peptide derived from human LIPJ AA range: 228-278 |
| Specificity | This antibody detects endogenous levels of human LIPJ |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C . Avoid repeated freeze-thaw cycles. |
| Protein Name | Lipase member J (EC 3.1.1.-) (Lipase-like abhydrolase domain-containing protein 1) |
| Gene Name | LIPJ LIPL1 |
| Cellular localization | |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | 142910 |
| Human Swiss-Prot Number | Q5W064 |
| Alternative Names | |
| Background | |

