



U2AFL rabbit pAb

Cat No.:ES19533

For research use only

Overview

Product Name	U2AFL rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human U2AFL AA range: 131-181
Specificity	This antibody detects endogenous levels of human U2AFL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 (CCCH type zinc finger, RNA-binding motif and serine/arginine rich protein 1) (U2(RNU2) small nuclear RNA auxiliary
Gene Name	ZRSR1 U2AF1-RS1 U2AF1L1 U2AF1P U2AF1RS1 U2AFBPL
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q15695
Alternative Names	
Background	

