

ENL rabbit pAb

Cat No.:ES19534

For research use only

Overview

Product Name	ENL rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ENL AA range: 37-87
Specificity	This antibody detects endogenous levels of human ENL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein ENL (YEATS domain-containing protein 1)
Gene Name	MLLT1 ENL LTG19 YEATS1
Cellular localization	Nucleus.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	4298
Human Swiss-Prot Number	Q03111
Alternative Names	
Background	

