

1A1L1 rabbit pAb

Cat No.: ES19540

For research use only

Overview

Product Name 1A1L1 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human 1A1L1 AA

range: 242-292

Specificity This antibody detects endogenous levels of human

1A1L1

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name 1-aminocyclopropane-1-carboxylate synthase-like

protein 1 (ACC synthase-like protein 1)

Gene Name ACCS PHACS

Cellular localization

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/ml

Observed band

Human Gene ID 84680 Human Swiss-Prot Number Q96QU6

Alternative Names

Background

