



VA0E2 rabbit pAb

Cat No.:ES19542

For research use only

Overview

Product Name	VA0E2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human VA0E2 AA range: 26-75
Specificity	This antibody detects endogenous levels of human VA0E2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	V-type proton ATPase subunit e 2 (V-ATPase subunit e 2) (Lysosomal 9 kDa H(+)-transporting ATPase V0 subunit e2) (Vacuolar proton pump subunit e 2)
Gene Name	ATP6V0E2 ATP6V0E2L C7orf32
Cellular localization	Membrane ; Multi-pass membrane protein . Cytoplasmic vesicle, clathrin-coated vesicle membrane ; Multi-pass membrane protein . Cytoplasmic vesicle, secretory vesicle, synaptic vesicle membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	155066
Human Swiss-Prot Number	Q8NHE4
Alternative Names	
Background	

