

# MCMD2 rabbit pAb

Cat No.:ES19543

For research use only

## Overview

Product Name	MCMD2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MCMD2 AA range: 488-538
Specificity	This antibody detects endogenous levels of human MCMD2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	MCM domain-containing protein 2
Gene Name	MCMD2 C8orf45
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	157777
Human Swiss-Prot Number	Q4G0Z9
Alternative Names	
Background	

