



SSPN rabbit pAb

Cat No.:ES19545

For research use only

Overview

Product Name	SSPN rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SSPN AA range: 154-204
Specificity	This antibody detects endogenous levels of human SSPN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Sarcospan (K-ras oncogene-associated protein) (Kirsten-ras-associated protein)
Gene Name	SSPN KRAG
Cellular localization	Cell membrane ; Multi-pass membrane protein . Cell membrane, sarcolemma . Cell junction, synapse, postsynaptic cell membrane ; Multi-pass membrane protein . Also found in myotendinous junctions and in the postsynaptic membrane of neuromuscular junctions. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	8082
Human Swiss-Prot Number	Q14714
Alternative Names	
Background	

