



DOCK6 rabbit pAb

Cat No.:ES19546

For research use only

Overview

Product Name	DOCK6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DOCK6 AA range: 1460-1510
Specificity	This antibody detects endogenous levels of human DOCK6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Dedicator of cytokinesis protein 6
Gene Name	DOCK6 KIAA1395
Cellular localization	Cytoplasm . Cytoplasm, perinuclear region . Mainly located near the cell surface.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	57572
Human Swiss-Prot Number	Q96HP0
Alternative Names	
Background	

