



IZUM3 rabbit pAb

Cat No.:ES19547

For research use only

Overview

Product Name	IZUM3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human IZUM3 AA range: 32-82
Specificity	This antibody detects endogenous levels of human IZUM3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Izumo sperm-egg fusion protein 3
Gene Name	IZUMO3 C9orf134
Cellular localization	Cell membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q5VZ72
Alternative Names	
Background	

