

# PARVG rabbit pAb

Cat No.:ES19548

For research use only

## Overview

Product Name	PARVG rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PARVG AA range: 113-163
Specificity	This antibody detects endogenous levels of human PARVG
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Gamma-parvin
Gene Name	PARVG
Cellular localization	Cell junction, focal adhesion. Cell membrane ; Peripheral membrane protein ; Cytoplasmic side . Cytoplasm, cytoskeleton . Constituent of focal adhesions. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	64098
Human Swiss-Prot Number	Q9HBI0
Alternative Names	
Background	

