



TLE6 rabbit pAb

Cat No.:ES19549

For research use only

Overview

Product Name	TLE6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TLE6 AA range: 144-194
Specificity	This antibody detects endogenous levels of human TLE6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Transducin-like enhancer protein 6
Gene Name	TLE6
Cellular localization	Nucleus .; [Isoform 1]: Cytoplasm .; [Isoform 2]: Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79816
Human Swiss-Prot Number	Q9H808
Alternative Names	
Background	

