



KDM4E rabbit pAb

Cat No.:ES19550

For research use only

Overview

Product Name	KDM4E rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KDM4E AA range: 338-388
Specificity	This antibody detects endogenous levels of human KDM4E
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Lysine-specific demethylase 4E (EC 1.14.11.-) (KDM4D-like protein) (Lysine-specific demethylase 4D-like)
Gene Name	KDM4E KDM4DL
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	390245
Human Swiss-Prot Number	B2RXH2
Alternative Names	
Background	

