

ATLA1 rabbit pAb

Cat No.:ES19558

For research use only

Overview

Product Name	ATLA1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ATLA1 AA range: 460-510
Specificity	This antibody detects endogenous levels of human ATLA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Atlastin-1 (EC 3.6.5.-) (Brain-specific GTP-binding protein) (GTP-binding protein 3) (GBP-3) (hGBP3) (Guanine nucleotide-binding protein 3) (Spastic paraplegia 3 protein A)
Gene Name	ATL1 GBP3 SPG3A
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Golgi apparatus membrane ; Multi-pass membrane protein . Cell projection, axon . Localizes to endoplasmic reticulum tubular network (PubMed:27619977). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	51062
Human Swiss-Prot Number	Q8WXF7
Alternative Names	
Background	

