

CHDH rabbit pAb

Cat No.:ES19559

For research use only

Overview

Product Name	CHDH rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CHDH
Specificity	This antibody detects endogenous levels of human CHDH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	Choline dehydrogenase, mitochondrial (CDH) (CHD) (EC 1.1.99.1)
Gene Name	CHDH
Cellular localization	Mitochondrion inner membrane .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55349
Human Swiss-Prot Number	Q8NE62
Alternative Names	
Background	

