



SYAM rabbit pAb

Cat No.:ES19560

For research use only

Overview

Product Name	SYAM rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SYAM AA range: 495-545
Specificity	This antibody detects endogenous levels of human SYAM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Alanine--tRNA ligase, mitochondrial (EC 6.1.1.7) (Alanyl-tRNA synthetase) (AlaRS)
Gene Name	AARS2 AARSL KIAA1270
Cellular localization	Mitochondrion .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	57505
Human Swiss-Prot Number	Q5JTZ9
Alternative Names	
Background	

