

MKLN1 rabbit pAb

Cat No.:ES19561

For research use only

Overview

Product Name	MKLN1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MKLN1 AA range: 61-111
Specificity	This antibody detects endogenous levels of human MKLN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Muskelin
Gene Name	MKLN1
Cellular localization	Cytoplasm . Cytoplasm, cytosol . Nucleus, nucleoplasm . Cell projection, ruffle . Cytoplasm, cell cortex . Cell junction, synapse . Cell junction, synapse, postsynapse . Colocalizes with GABRA1 at synapses and in postsynaptic regions. Colocalizes with actin fibers in the cell cortex. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	4289
Human Swiss-Prot Number	Q9UL63
Alternative Names	
Background	

