



RRMJ3 rabbit pAb

Cat No.:ES19564

For research use only

Overview

Product Name	RRMJ3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RRMJ3 AA range: 628-678
Specificity	This antibody detects endogenous levels of human RRMJ3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	pre-rRNA processing protein FTSJ3 (EC 2.1.1.-) (Protein ftsJ homolog 3) (Putative rRNA methyltransferase 3) (rRNA (uridine-2'-O-)-methyltransferase 3)
Gene Name	FTSJ3 SB92
Cellular localization	Nucleus, nucleolus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	117246
Human Swiss-Prot Number	Q8IY81
Alternative Names	
Background	

