



SO1B7 rabbit pAb

Cat No.:ES19566

For research use only

Overview

Product Name	SO1B7 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SO1B7 AA range: 270-320
Specificity	This antibody detects endogenous levels of human SO1B7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Putative solute carrier organic anion transporter family member 1B7 (Liver-specific organic anion transporter 3) (LST-3)
Gene Name	SLCO1B7 LST3 LST3TM12
Cellular localization	Cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	338821
Human Swiss-Prot Number	G3V0H7
Alternative Names	
Background	

